

## CD326 S-pluriBead® anti-ms

**Order No.:**29-32600-10



### 概要

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

### Additional Information

|        |   |
|--------|---|
| 大きさ    | 2 ml  |
| しゅるい   | マウス   |
| 細胞型    | 循環腫瘍細胞, CTC (循環腫瘍細胞, 上皮細胞)  |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔     | >90%  |
| 保管状態   | 4-8 °C  |
| きせいぶん  | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.