

universal S-pluriBead® anti-hu
Order No.:31-RBaHU-10

概要

pluriBead® universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

Additional Information

| | |
|--------|---|
| 大きさ | 2 ml |
| しゅるい | 人間 |
| 細胞型 | 何れかの |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔 | >90% |
| 保管状態 | 4-8 °C |
| きせいぶん | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.