

CD47 S-pluriBead® anti-hu

Order No.:19-04700-10



概要

CD47 pluriBead® suspension for the positive cell isolation of CD47+ target cells from any human sample material. (2ml)

Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	血小板, 白血球, 赤血球, 赤血球
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.