

CD4 S-pluriBead® anti-hu

Order No.:19-00400-10



概要

CD4 pluriBead® suspension for the positive cell isolation of CD4+ target cells from any human sample material. (2ml)

Additional Information

| | |
|--------|---|
| 大きさ | 2 ml |
| しゅるい | 人間 |
| 細胞型 | 単球 (サブセット), ヘルパーT細胞 |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔 | >90% |
| 保管状態 | 4-8 °C |
| きせいぶん | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.