

## CD7 S-pluriBead® anti-hu

**Order No.:**19-00700-10



### 概要

CD7 pluriBead® suspension for the positive cell isolation of CD7+ target cells from any human sample material. (2ml)

### Additional Information

|        |   |
|--------|---|
| 大きさ    | 2 ml  |
| しゅるい   | 人間  |
| 細胞型    | NK細胞, T細胞, 胸腺細胞   |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔     | >90%  |
| 保管状態   | 4-8 °C  |
| きせいぶん  | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.