


10% BSA stock solution in 1x PBS (500 ml)

Order No.:60-00020-12



概要

10% BSA stock solution in 1x PBS can be used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

Additional Information

大きさ	500 ml
納期	1-2
興味領域	細胞生物学
不妊症	無菌状態
フオーンレンジ	pH 7.4
うんぱんじょうたい	室温
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.