


10x PBS Solution



Order No.	Size	Sterility
60-00010-10	100 ml	無菌状態
60-00010-11	500 ml	無菌状態

概要

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without**:

- Calcium
- Magnesium
- Phenol Red
- Protease free

Components

- NaCl (1368,9 mMol/L)
- Na₂HPO₄ • 2H₂O (100mMol/L)
- KH₂PO₄ (19,8 mMol/L)
- KCl (26,8 mMol/L)

Additional Information

納期	1-2
興味領域	細胞生物学
不妊症	無菌状態
フォーレンジ	pH 6.8
うんぱんじょうたい	室温
保管状態	室温

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.