

EDTA Solution (0.2M in PBS)



Order No.	Size	Sterility
60-00030-10	10 ml	無菌状態
60-00030-11	50 ml	無菌状態
60-00030-12	100 ml	無菌状態

概要

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

興味領域	細胞生物学
不妊症	無菌状態
フォーレンジ	pH 7.4
うんぱんじょうたい	室温
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.