

pluriSpin® Human PBMC Depletion (2 ml)

Order No.:19-08002-31



概要

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	非PBMC細胞
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.