

## pluriSpin® Human CD19<sup>+</sup> Cell Depletion



Order No.	Size
19-06002-31	2 ml
19-06002-32	10 ml

### 概要

pluriSpin® Human CD19<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched b cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels B cells, while the isolated human b cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified b cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

しゅるい	人間
細胞型	B 細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.