

## pluriSpin® Human CD4<sup>+</sup> Cell Depletion





| Order No.   | Size  |
|-------------|-------|
| 19-03002-31 | 2 ml  |
| 19-03002-32 | 10 ml |

### 概要

pluriSpin® Human CD4<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched helper t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Helper T cells, while the isolated human helper t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified helper t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|        |   |
|--------|---|
| しゅるい   | 人間  |
| 細胞型    | ヘルパーT細胞   |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔     | >90%  |
| 保管状態   | 室温  |
| きせいぶん  | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.