

pluriSpin® Human CD45⁺ Cell Depletion



Order No.	Size
19-09002-31	2 ml
19-09002-32	10 ml

概要

pluriSpin® Human CD45⁺ Cell Depletion is designed to isolate pure, viable and untouched CD45⁺-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45⁺ cells, while the isolated population free of CD45⁺ are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45⁺ cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

しゅるい	人間
細胞型	CD45-細胞
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.