

pluriSpin® Human Monocyte Depletion


Order No.
Size

19-01002-31

2 ml

19-01002-32

10 ml

概要

pluriSpin® Human Monocyte Depletion is designed to isolate pure, viable and untouched monocytes-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels monocytes, while the isolated human monocytes-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

しゅるい	人間
細胞型	単球
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.