

pluriSpin® Human PBMC Depletion


Order No.
Size

19-08002-31

2 ml

19-08002-32

10 ml

概要

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

しゅるい	人間
細胞型	非PBMC細胞
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.