

## ÜberStrainer 1 µm



**Order No.:**43-70001-03

### 概要

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

### Additional Information

メッシュサイズ	1 µm
大きさ	10 pcs
ファブリック素材	PET (ポリエチルテレフタレート)
ハウジング材	PP (ポリプロピレン)
Mesh Type	woven
不妊症	無菌状態
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

