

Bundle pluriStrainer® Mini 5 µm and Flow Cytometry Tubes (25 pcs.)

Order No.:43-10005-46



概要

Bundle includes sterile cell strainers - pluriStrainer Mini with 5 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 5 µm polyethylenterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 5 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

アプリケーション

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

Additional Information

納期	1-2
メッシュサイズ	5 µm
包装	25 pcs. (25 pcs. / Bag)
大きさ	25 pcs.
色	黄色い
フローコントロール	いいえ
Mesh Fixation	Injected in housing for strong hold
ファブリック素材	PET (ポリエチルテレフタレート)

pluriSelect Life Science

Deutscher Platz 5a, 04103 Leipzig, DE
 Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99
 support@pluriselect.com
<https://www.pluriselect.com/>

ハウジング材	PP (ポリプロピレン)
嵌まる	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
扱い	Gripping surface enables easy and aseptic handling.
不妊症	無菌状態
に匹敵する	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
安定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
えんしんかのう	はい
組織解離	はい
しょうどくかのう	はい
うんぱんじょうたい	室温
保管状態	室温
ほうていじょうほう	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.