

## Bundle pluriStrainer® Mini 20 µm and Flow Cytometry Tubes (500 pcs.)

**Order No.:**43-10020-76



### 概要

Bundle includes sterile pluriStrainers Mini with 20 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 20 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 20 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

### アプリケーション

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marcol-sm-12, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry(FACS™)

### Additional Information

納期	1-2
メッシュサイズ	20 µm
包装	500 pcs. (25 pcs. / Bag)
大きさ	500 pcs.
色	青々
フローコントロール	いいえ
Mesh Fixation	Injected in housing for strong hold
ファブリック素材	PET (ポリエチルテレフタレート)

ハウジング材	PP (ポリプロピレン)
嵌まる	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
扱い	Gripping surface enables easy and aseptic handling.
不妊症	無菌状態
に匹敵する	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
安定性	RLT® Buffer, TRIzol®, isopropanol, organic solvents
えんしんかのう	はい
組織解離	はい
しょうどくかのう	はい
うんぱんじょうたい	室温
保管状態	室温
ほうていじょうほう	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

## Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.