

## ÜberStrainer Set 1

**Order No.:**43-70000-97



### 概要

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 20 µm, 60 µm, 100 µm, 300 µm, 500 µm.

### Additional Information

メッシュサイズ	Multiple
大きさ	10 pcs
ファブリック素材	PET (ポリエチルテレフタレート)
ハウジング材	PP (ポリプロピレン)
Mesh Type	woven
不妊症	無菌状態
うんぱんじょうたい	室温
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures.

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.