

ÜberStrainer 400 µm



Order No.:43-70400-03

概要

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

Additional Information

メッシュサイズ	400 µm
大きさ	10 pcs
ファブリック素材	PET (ポリエチレンテレフタレート)
ハウジング材	PP (ポリプロピレン)
Mesh Type	woven
不妊症	無菌状態
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

