

## Re-Strainer (In-Line-Strainer) 1 µm - Sterile

**Order No.:**43-75001-40



### 概要

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### アプリケーション

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Additional Information

納期	1-2
メッシュサイズ	1 µm
変化	20 pcs, sterile
包装	4 devices per bag

大きさ	20 pcs.
色	あか
フローコントロール	はい
Mesh Fixation	Injected in housing for strong hold
ファブリック素材	PET (ポリエチルテレフタレート)
ハウジング材	LD-PE (低密度ポリエチレン)
不妊症	無菌状態
安定性	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
うんぱんじょうたい	Room Temperature
保管状態	室温
きせいぶん	For research use only
ほうていじょうほう	<b>The product is for research and development only, not for diagnostic or therapeutic use</b>

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.