

Re-Strainer (In-Line-Strainer) 100 µm - Sterile

Order No.:43-75100-40



概要

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

アプリケーション

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

Additional Information

納期	1-2
メッシュサイズ	100 µm
包装	4 devices per bag
大きさ	20 pcs.

色	黄色い
フローコントロール	はい
Mesh Fixation	Injected in housing for strong hold
ファブリック素材	PET (ポリエチルテレフタレート)
ハウジング材	LD-PE (低密度ポリエチレン)
不妊症	無菌状態
安定性	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
うんぱんじょうたい	Room Temperature
保管状態	室温
きせいぶん	For research use only
ほうていじょうほう	The product is for research and development only, not for diagnostic or therapeutic use

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.