

pluriSpin® Human CD4⁺ Cell Depletion (10 ml)

Order No.:19-03002-32



概要

pluriSpin® Human CD4⁺ Cell Depletion is designed to isolate pure, viable and untouched helper t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Helper T cells, while the isolated human helper t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified helper t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大きさ	10 ml
しゅるい	人間
細胞型	ヘルパーT細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

