

pluriSpin® Human CD8⁺ T Cell Enrichment (10 ml)

Order No.:19-04001-32



概要

pluriSpin® Human CD8⁺ T Cell Enrichment is designed to isolate pure, viable and untouched cytotoxic T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, helper T cells, B cells, NK cells and granulocytes, while the negatively isolated human cytotoxic T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic T cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大きさ	10 ml
しゅるい	人間
細胞型	細胞障害性T細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

