

pluriSpin® Human CD8⁺ Cell Depletion (2 ml)

Order No.:19-04002-31



概要

pluriSpin® Human CD8⁺ Cell Depletion is designed to isolate pure, viable and untouched cytotoxic t cell-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Cytotoxic T cells, while the isolated human cytotoxic t cell-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic t cell-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大きさ	2 ml
しゅるい	人間
細胞型	細胞障害性T細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

