

## pluriSpin® Human T Cell Depletion


**Order No.**
**Size**

19-02002-31

2 ml

19-02002-32

10 ml

**概要**

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

**Additional Information**

納期	1-2
しゅるい	人間
細胞型	T細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	4-8 °C
きせいぶん	For research use only. Not for use in diagnostic procedures,

**Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.

