

## pluriSpin® Human Monocyte Enrichment (2 ml)

**Order No.:**19-01001-31



### 概要

pluriSpin® Human Monocyte Enrichment is designed to isolate pure, viable and untouched monocytes from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, B cells, NK cells and granulocytes, while the negatively isolated human monocytes are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|        |   |
|--------|---|
| 大きさ    | 2 ml  |
| しゅるい   | 人間  |
| 細胞型    | 単球  |
| サンプル材料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 純潔     | >90%  |
| 保管状態   | 4-8 °C  |
| きせいぶん  | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

