

## 10x Wash Buffer (500 ml)

Order No.:60-00080-11



### 상세 설명

10x Wash Buffer (containing PBS/EDTA/BSA) - for washing of any sample material containing cells

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 500 ml   |
| Storage Condition    | 4-8 °C   |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.