

S-pluriBead® Maxi Reagent Kit

Order No.:70-50010-12



상세 설명

Maxi S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 S-pluriStrainer, 10 Funnel, 10 connector rings

Additional Information

| | |
|----------------------|---|
| Size | 10 |
| Delivery Time (days) | 1-2 |
| Species | Human |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | Room Temperature; Buffer 4-8 °C; Buffer D -20°C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.