

## CD47 M-pluriBead® anti-hu

**Order No.:**19-04700-20



### 상세 설명

CD47 pluriBead® suspension for the positive cell isolation of CD47+ target cells from any human sample material. (2ml)

### Additional Information

Size	2 ml
Species	Human
Cell Types	Platelets, Leukocytes, Erythrocytes, RBC
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.