

## CD19 M-pluriBead® anti-hu



**Order No.:**19-01900-20

### 상세 설명

CD19 pluriBead® suspension for the positive cell isolation of CD19+ target cells from any human sample material. (2ml)

### Additional Information

Size	2 ml
Species	Human
Cell Types	B cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.