

PLT Spin Medium (500 ml)

Order No.: 60-00094-12



상세 설명

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Additional Information

Size	500 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.