

## CD43 M-pluriBead® anti-hu

**Order No.:**19-04300-20



### 상세 설명

CD43 pluriBead® suspension for the positive cell isolation of CD43+ target cells from any human sample material. (2ml)

### Additional Information

Size	2 ml
Species	Human
Cell Types	Monocytes, T cells, Granulocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Shipping Condition	Room Temperature
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.