

CD326 S-pluriBead® anti-ms

Order No.:29-32600-10



상세 설명

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Circulating Tumor Cells, CTC, Epithelial cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.