

pluriSpin® Human PBMC Depletion (2 ml)

Order No.:19-08002-31



상세 설명

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Human |
| Cell Types | Non-PBMC cells |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.