

## pluriSpin® Human B Cell Enrichment





Order No.	Size
19-06001-31	2 ml
19-06001-32	10 ml

### 상세 설명

pluriSpin® Human B Cell Enrichment is designed to isolate pure, viable and untouched B cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, NK cells and granulocytes, while the negatively isolated human B cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified B cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Species	Human
Cell Types	B cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>80%
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.