

## pluriSpin® Human CD45<sup>+</sup> Cell Depletion





| Order No.   | Size  |
|-------------|-------|
| 19-09002-31 | 2 ml  |
| 19-09002-32 | 10 ml |

### 상세 설명

pluriSpin® Human CD45<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched CD45<sup>+</sup>-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45<sup>+</sup> cells, while the isolated population free of CD45<sup>+</sup> are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45<sup>+</sup> cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|                      |   |
|----------------------|---|
| Species              | Human   |
| Cell Types           | CD45 <sup>-</sup> Cell  |
| Storage Condition    | Room Temperature  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.