

pluriSpin® Human NK Cell Depletion



Order No.	Size
19-05002-31	2 ml
19-05002-32	10 ml

상세 설명

pluriSpin® Human NK Cell Depletion is designed to isolate pure, viable and untouched nk cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels NK cells, while the isolated human nk cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified nk cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Species	Human
Cell Types	NK cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.