

pluriSpin® Human PBMC Depletion (10 ml)

Order No.:19-08002-32



상세 설명

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

Size	10 ml
Species	Human
Cell Types	Non-PBMC cells
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.