

## Syringe Strainer 500 µm (12 pcs., non-sterile)



**Order No.:** 43-71500-50

### 상세 설명

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Mesh Size	500 µm
Variation	12 pcs, non-sterile
Size	12 pcs.
Flow Control	네]
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Non-Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.