

Monocytes Spin Medium


Order No.
Size

| | |
|-------------|--------|
| 60-00095-10 | 100 ml |
| 60-00095-11 | 250 ml |
| 60-00095-12 | 500 ml |

상세 설명

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

Additional Information

| | |
|----------------------|---|
| Delivery Time (days) | 1-2 |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Storage Condition | Room temperature, protect from light |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range | 6.8 – 7.4 |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.