

High Density Spin Medium 500 ml (1.1 g/ml)

Order No.:60-01103-12



상세 설명

Density Diluent Medium is a ready to use, sterile medium for the preparation of different density working solutions for the enrichment of cells, bacteria and other particles. The working solutions can be prepared by the dilution of High Density Spin Medium with Density Diluent Medium.

Additional Information

Size	500 ml
Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.