

CD5 S-pluriBead® anti-hu

Order No.:19-00500-10



Description

CD5 pluriBead® suspension for the positive cell isolation of CD5+ target cells from any human sample material. (2ml)

Additional Information

| | |
|----------------------|-------------------------------------------------------------------------------|
| Size | 2 ml |
| Species | Human |
| Cell Types | T cells, B cells (subset), Thymocytes |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.