

## pluriSpin® Human PBMC Depletion (2 ml)

**Order No.:**19-08002-31



### Description

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

### Additional Information

Size	2 ml
Species	Human
Cell Types	Non-PBMC cells
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.