

## Leuko Spin Medium (100 ml)

Order No.:60-00091-10



### Description

Leuko Spin is a sterile, ready-to-use density medium for the isolation of all Leukocytes from fresh human peripheral blood or buffy coat with a high yield. Mononuclear and polymorphonuclear (multi-lobed nuclei) cells can be enriched by a single-step density gradient centrifugation with Leuko Spin Medium. Through combination with PBMC Spin in a double gradient centrifugation it is possible to isolate and separate granulocytes and PBMCs in one step with high yield and purity. The enriched cells can be used for a wide range of downstream applications such as MACS (magnetic activated cell sorting, a trademark of Miltenyi Biotec GmbH) or stem cell research.

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 100 ml   |
| Delivery Time (days) | 1-2  |
| Storage Condition    | Room temperature, protect from light                         |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range             | 6.8 - 7.4  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.