

## Syringe Strainer 60 µm (12 pcs., non-sterile)

**Order No.:**43-71060-50



### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

|                      |                                                              |
|----------------------|--------------------------------------------------------------|
| Mesh Size            | 60 µm                                                        |
| Variation            | 12 pcs, non-sterile                                          |
| Size                 | 12 pcs.                                                      |
| Flow Control         | Yes                                                          |
| Housing Material     | LD-PE (Low Density Polyethylen)                              |
| Mesh Type            | woven                                                        |
| Sterility            | Non-Sterile                                                  |
| Shipping Condition   | Room Temperature                                             |
| Storage Condition    | Room Temperature                                             |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.