

## S-pluriBead® Maxi Reagent Kit

**Order No.:**70-50010-12



### Description

Maxi S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 S-pluriStrainer, 10 Funnel, 10 connector rings

### Additional Information

|                      |   |
|----------------------|---|
| Size                 | 10  |
| Delivery Time (days) | 1-2   |
| Species              | Human   |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%  |
| Storage Condition    | Room Temperature; Buffer 4-8 °C; Buffer D -20°C                               |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.