

## Leuko Spin Medium (500 ml)

**Order No.:** 60-00091-12



### Description

Leuko Spin is a sterile, ready-to-use density medium for the isolation of all Leukocytes from fresh human peripheral blood or buffy coat with a high yield. Mononuclear and polymorphonuclear (multi-lobed nuclei) cells can be enriched by a single-step density gradient centrifugation with Leuko Spin Medium. Through combination with PBMC Spin in a double gradient centrifugation it is possible to isolate and separate granulocytes and PBMCs in one step with high yield and purity. The enriched cells can be used for a wide range of downstream applications such as MACS (magnetic activated cell sorting, a trademark of Miltenyi Biotec GmbH) or stem cell research.

### Additional Information

Size	500 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.