

## 10x PBS Solution



Order No.	Size	Sterility
60-00010-10	100 ml	Sterile
60-00010-11	500 ml	Sterile

### Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without**:

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (1368,9 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (100mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (19,8 mMol/L)
- KCl (26,8 mMol/L)

### Additional Information

Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 6.8
Shipping Condition	Room Temperature
Storage Condition	Room Temperature

pluriSelect Life Science

Deutscher Platz 5a, 04103 Leipzig, DE

Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99

support@pluriselect.com

<https://www.pluriselect.com/>

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.