

Streck - Cell-Free DNA BCT® (10 ml, CE-IVD)



Order No.
Size

218996

6 tubes

218997

100 tubes

230244

1000 tubes

Description

A blood collection tube which stabilizes nucleated blood cells. The unique preservative prevents the release of genomic DNA, allowing isolation of high-quality cell-free DNA.

Convenient

- Reduces the need for immediate plasma preparation and CTC processing
- Allows for isolation of high-quality cell-free DNA for a wide range of downstream applications
- Allows for convenient sample collection, transport and storage

Stable

- cfDNA and gDNA are stable for up to 14 days at 6 °C to 37 °C
- CTCs are stable for up to 7 days at 15 °C to 30 °C
- Unfilled tubes have 24-month shelf life from the date of manufacture

Cell-Free DNA BCT is the trusted global leader in cell-stabilization

Additional Information

Prefilled	Unfilled
Delivery Time (days)	1-2
Storage Condition	Room Temperature

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

