

pluriSpin® Human T Cell Depletion



| Order No. | Size |
|-------------|-------|
| 19-02002-31 | 2 ml |
| 19-02002-32 | 10 ml |

Description

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Delivery Time (days) | 1-2 |
| Species | Human |
| Cell Types | T cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

